

primed



for the Perfect *Swipe!*

primedENDO.com

Manufacturer Of Top Quality Cleaning Brushes



Medical Device Reprocessing Starts with PriMed Brushes!



• Rigid Instruments Brushes

MODEL NUMBER	DESCRIPTION
88816	Double Ended Instrument Cleaning Brush
88800	Large Instrument Cleaning Brush 7"
88303	Instrument Cleaning Brush 7"
88755	Acetabular Reamer Brush
88905	Baron Tube Cleaning Brush – 5FR (1.6mm)
88906	Frazier Tube Cleaning Brush – 6FR (2mm)
88907	Frazier Tube Cleaning Brush – 7FR (2.3mm)
88908	Frazier Tube Cleaning Brush – 8FR (2.7mm)
88909	Frazier Tube Cleaning Brush – 9FR (3mm)
88910	Frazier Tube Cleaning Brush – 10FR (3.3mm)
88911	Frazier Tube Cleaning Brush – 11FR (3.7mm)
88912	Frazier Tube Cleaning Brush – 12FR (4mm)
88539	Twisted Wire Cleaning Brush – 5mm 16"
88540	Twisted Wire Cleaning Brush – 7mm 16"
88541	Twisted Wire Cleaning Brush – 10mm 16"
88542	Twisted Wire Cleaning Brush – 2mm x 12"
88543	Twisted Wire Cleaning Brush – 3mm x 12"
88544	Twisted Wire Cleaning Brush – 4mm x 12"
88545	Twisted Wire Cleaning Brush – 5mm x 12"
88546	Twisted Wire Cleaning Brush – 6mm x 12"
88550	Twisted Wire Cleaning Brush – 2.4mm x 24"
88551	Twisted Wire Cleaning Brush – 5mm x 24"
88552	Twisted Wire Cleaning Brush – 7mm x 24"
88553	Twisted Wire Cleaning Brush – 10mm x 24"
85001	Twisted Wire Cleaning Brush – 5mm x 24" Stiff Shaft Blue bristles with protective bullet tip
85002	Twisted Wire Cleaning Brush – 5mm x 24" Stiff Shaft Blue bristles
85003	Twisted Wire Cleaning Brush – 8mm x 24" Stiff Shaft Blue bristles with protective bullet tip



• Brush Hanger Available - Model Number BH20

• Single Use Flexible Endoscopes Cleaning Brushes

SINGLE USE CHANNEL CLEANING BRUSHES	
MODEL NUMBER	DESCRIPTION
9190	Minimum channel size 2 mm. Length 240 cm, 4-6.3 mm tapered brush head.
9192 Triple Head	Minimum channel size 2 mm. Length 240 cm, 4-6.3 mm tapered brush heads. Push & Pull.
9194 Double Ended	Minimum channel size 2 mm. Length 240 cm, 4-6.3 mm tapered brush heads. Push & Pull.
9123 Double Ended	Minimum channel size 2 mm. Length 120 cm, 3 mm brush heads. Push & Pull.
9115	Minimum channel size 2 mm. Length 115 cm, 3-4 mm tapered brush head.
9025 Micro Brush	Minimum channel size 1.2 mm. Length 120 cm, 2.5 mm brush head.
9020 Micro Brush	Minimum channel size 1.2 mm. Length 120 cm, 2 mm brush head.
9010 Micro Brush	Minimum channel size 1 mm. Length 110 cm, 2.5 mm brush head.
9103 Micro Brush	Minimum channel size 1 mm. Length 220cm, 2.5 mm brush head.

SINGLE USE SUCTION WELL AND VALVE CLEANING BRUSHES	
MODEL NUMBER	DESCRIPTION
9990 Suction Well	Working Length 5.5 cm, 11 mm brush head.
9992 Double Ended Combo	Working Length 5.5 cm & 5.5 cm, 4.3 & 11 mm brush heads.
99PTX2 Double Ended Combo	Working Length 5.5 cm & 5.5 cm, 5.8 & 10-12 mm tapered brush heads.

SINGLE USE ENDOSCOPE CLEANING BRUSH KITS	
MODEL NUMBER	DESCRIPTION
9199	Channel & Suction Well Cleaning Brush Kit: (9190 + 9990) - 2 pieces
9155	Channel & Suction Well/Valve Cleaning Brush Kit: (9115 + 9992) - 2 pieces
9142	Channel & Suction Well/Valve Cleaning Brush Kit: (9194 + 9992) - 2 pieces
9140	Channel & Suction Well Cleaning Brush Kit: (9194 + 9990) - 2 pieces
9122	Channel & Suction Well/Valve Cleaning Brush Kit: (9192 + 9992) - 2 pieces

- Other sizes, custom configurations and combinations are available.
- Please contact us with your requirements.



ORDERING INFORMATION

Email: sales@primedENDO.com

Phone Toll Free: **877-565-0565**

Fax Toll Free: **877-565-0566**

Manufactured by:
PriMed Instruments Inc.
1080 Tristar Drive, Unit #14
Mississauga, ON L5T 1P1
Canada

905-565-0565

905-565-0566